



AM7 Electric Motor AI Motor Series

LED

Rotation Speed: 2000-40000 rpm

Water Spray: Internal water

Voltage input: DC 24V

Micro-motor Size: $\Phi 22 \times 76.7$ mm

Micro-motor Material: Stainless

Main unit Size : 120X145X50mm

Air Prussure : (2.5 bar ~ 5 bar), air flow > 1.5L / min

Connect : M4 / B2 / K / N

Water Prussure : (2 bar ~ 5 bar), water flow > 50ml / min

• With LED(white)

Product List

	Main unit (1 PCS)		Micro-motor (1 PCS)
	O-ring (4 PCS)		Disinfecting Plug (1 PCS)
	Disinfecting Cap (1 PCS)		AC adapter (1 PCS)
	Input Power Cord (1 PCS)		Manual (1 PCS)



Packaging Display & Size:

Material: Paper+Sponge

Size:..... 275X205X115mm

Weight: 140g $\pm 10\%$