

AM8 AI Motor Series

LED

Rotation Speed: 1000-40000 rpm

Dual Mode : Surgery&Restoration

Water Spray: Internal water/External Water

Micro-motor Size: $\Phi 22 \times 76.7 \text{ mm}$

Micro-motor Material: Stainless

Main unit Size : 175X190X230mm

Air Prussure : (2.5 bar ~ 5 bar), air flow > 1.5L / min

Water Prussure : (2 bar ~ 5 bar), water flow > 50ml / min

• With LED(white)

Product List

	Main unit (1 PCS)		Micro-motor (1 PCS)
	SG1-4.2L (1 PCS)		S1-5L (1 PCS)
	Input Power Cord (1 PCS)		AC adapter (1 PCS)
	Footstep (1 PCS)		Manual (1 PCS)
	Disinfecting Cap (1 PCS)		Disinfecting Plug (1 PCS)
○○○○○		O-ring (4 PCS)	

