

S670

DENTAL UNIT

SIMPLE AND COMFORTABLE
PRACTICAL AND RELIABLE



S670

DENTAL UNIT



ENGINEERING ELECTRIC CHAIR

Comfortable chair cushion design

The seat design is based on the human body shape and provides perfect comfort and support, especially for patients who have been treated for a long time.



S670

DENTAL UNIT



- **LED Lamp**
The handle is removable for cleaning, sterilization



- **Sensor on/off** (white light, yellow light and mixed 3 kinds of light source mode manual conversion) Reduces eye fatigue and precise colorimetry



- **Assistant holder combination (off delay)**
Aluminum handle can be removed for cleaning and disinfection




- **Assistant holder silicone pad**
Scratch resistant, easy to clean and disinfectable.




- **Aluminum instrument plate (Interlock)**
The handle is removable for disinfection and the hanger is rotatable.



- **Instrument tray silicone pad**
Scratch resistant, easy to clean and disinfectable.




- **External water storage bottle**
One-key conversion, one-key switch, convenient and fast.



- **Removable strength and weakness suction filtration device** easy to clean

S670

DENTAL UNIT



- **Double joint headrest**
The fold-out design makes the patient more comfortable and at ease during treatment



- **Silent electric chair (Security Protection)**
The lowest chair position is 380mm design, which is convenient for children and elderly patients



- **Rotatable handrail design**
Aluminium alloy/engineering plastic face cover



- **Emergency stop switch**
One-click stop, safe and convenient



- **Ceramic spittoon, automatic water supply system**
Removable/rotating, rinsing mouth ration water.



- **Floor-standing box**
Installed on the ground, high stability.



- **Multifunctional foot pedal**
cup/water supply/oral lights/ spitting position/handpiece working water Control up and down/supine/flushing



- **Water/air/electricity switchgear (Anti-suction device)**
One-click on/off convenient and fast.

S1 LED LAMP	☐
S3 LED LAMP	☐
Disinfectior for water	☐
Built out floorbox	☐
Inductive water supply	☐
Glassspittoon	☐
Intraoral camera(with bracket))	☐
Central suction machine	☐
Down-mounted instrument tray	☐
Cart type instrument tray	☐
Teeth scaler machine(built-in)	☐
LED light curing machine(built-in)	☐
Fiber leather	☐
Imported/domestic high-speed handpiece	☐
Imported/domestic low-speed handpiece	☐
Imported/domestic electric motor(built-in)	☐